



Plastomed

Systems for your life...



SUCTION, BUBBLE & SILICON TUBES



Innovative plastic products and systems for medicine and technology

Product description

As one of the leading European manufacturers of suction hoses, we produce a comprehensive and flexible range of PVC suction hoses for single use applications.

This wide product range results from different possible combinations of:

- a wide variety of hose connectors (connectors and funnels)
- PVC materials, design as standard hose Shore A70 or softer hose in Shore A65
- diameter in various Charrière sizes (CH 25 and 30 as well as intermediate sizes on request)
- individual hose lengths, possible up to 600 cm
- sterile or non-sterile
- individually packed in peel packaging (sterile/non-sterile) or individually packed non-sterile in polyethylene bags

Therefore we can offer you a wide variety of suction hoses for your vacuum system adapted to your individual requirements and specifications.

All suction hoses are supplied individually packaged as a peel packaging in a sterile or non-sterile version.

Alternatively, we can also package non-sterile suction hoses individually in PE bags. You can then assemble our products into individual surgery sets, pack them and sterilise them.

Suction hoses must remain free of kinks to ensure their proper function.

The balance of flexibility and stability required for this is achieved in our suction hoses by an optimal diameter-wall thickness ratio.

All Plastomed suction hoses are made of PVC, are free of latex components and contain phthalate-free softeners. All suction hoses are also available as soft hoses (softer material).

Our hose connectors ensure a secure connection or disconnection of the suction hose and are suitable for connection to the collection containers of all vacuum systems. They can be connected to any end pieces such as suction catheters, drains, curettes and suction and suction cannulas and suction nozzles.

Notes when ordering

On the following pages we give you an overview of our standard product range.

In addition, we can also manufacture customised suction hoses on request. Please do not hesitate to contact us.

Beispiel Artikelzusammenstellung:

Art.-Nr.:	AS-XXX-CH TFS
AS:	Absaugschlauch
XXX:	Länge des Schlauches in cm
CH:	Charrière Größe des Schlauches
TF:	Schlauchverbindungsstücke, z.B. Universaltrichter (T) und Fingertip (F)
S:	Steril (ohne S = unsteril)

Contents

Tube attachments/ Variations		page
Universal connectors		4
Stepped conical connectors		4
Anti-kink conical connectors		4
Catheter connectors		4
Straight attachments		5
Fingertip		5
Stepped connector / Conus		5
Suction tubes		
Suction tubes with universal connector and fingertip		6
Suction tubes with universal connector		6
Suction tubes with fingertip		6
Suction tubes without attachments		7
Suction tubes with 2 universal connectors		7
Suction tubes with universal connector and loose fingertip		7
Suction tubes with 2 straight attachments		8
Suction tubes with anti-kink conical connector and catheter connector		8
Suction tubes with 2 anti-kink conical connectors		8
Suction tubes with catheter connector and fingertip		9
Suction tubes with 2 catheter connectors		9
Suction tubes with anti-kink conical connector and fingertip		9
Suction tubes with stepped conical connector and fingertip		10
Suction tubes with stepped conical connector and straight attachment		10
Suction tubes with 2 stepped conical connectors		10
Suction tubes with stepped conical connector and universal connector		11
Suction tubes with stepped conical connector		11
Suction tubes with stepped conical connector and anti-kink conical connector		11
Suction tubes with stepped conical connector and conus		12
Bubble tubes		12
Fingertip		12
Silicon tubes and sets		
Silicon tubes		14
Ordering information		
Order formula		15

TUBE ATTACHMENTS / VARIATIONS



UNIVERSAL CONNECTORS (T)

Varied application onto suction pumps or bags with small connections. Tubes with universal funnel find use in fields such as in home-care, intensive care stations, OP theatres, tracheotomy or in ENT.



STEPPED CONICAL CONNECTORS (Tg)

Varied application onto suction pumps or bags with small and large connections. The anti-kink feature ensures a premature collapse of the tube is prevented under loading, because the adapter bends. Tubes with stepped conical connector find use in areas such as in OP theatres.



ANTI-KINK CONICAL CONNECTORS (Tk)


Varied application onto suction pumps or bags with small and large connections. The anti-kink feature ensures a premature collapse of the tube is prevented under loading, because the adapter bends. Tubes with anti-kink conical connector find use in areas such as in OP theatres.



CATHETER CONNECTORS (Tb)

Varied application onto suction pumps or bags with small and large connections. Tubes with catheter adapter find use in areas such as in in home-care, intensive care stations, OP theatres, tracheotomy or ENT.

TUBE ATTACHMENTS / VARIATIONS



STRAIGHT ATTACHMENTS (L)

Varied application onto suction pumps or bags with small and large connections. Tubes with straight attachments find use in areas such as in OP theatres.



FINGERTIP (F)

Interrupter attachment (fingertip) purpose: regulation of the suction flow of fluids and OP residues.



STEPPED CONNECTOR / CONUS (K)

The connector attaches the suction tube as extension to the suction unit with the thorax drainage. Tubes with stepped connectors find application in OP theatres.



SCHLAUCHKONNEKTOR (R)

Schlauchverbindungsstück für Anschluss eines weiteren Absaug Schlauch, erhältlich für Absaugschläuche in den Grössen Charrière 30.
Auf Anfrage!


SUCTION TUBES

SUCTION TUBES WITH UNIVERSAL CONNECTOR (T) AND FINGERTIP (F)




REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TF	S	150	25	100
AS - 180 - 25 TF	S	180	25	100
AS - 200 - 25 TF	S	200	25	100
AS - 210 - 25 TF	S	210	25	100
AS - 240 - 25 TF	S	240	25	75
AS - 300 - 25 TF	S	300	25	60
AS - 350 - 25 TF	S	350	25	50
AS - 400 - 25 TF	S	400	25	50
AS - 150 - 30 TF	S	150	30	50
AS - 180 - 30 TF	S	180	30	50
AS - 200 - 30 TF	S	200	30	50
AS - 210 - 30 TF	S	210	30	50
AS - 240 - 30 TF	S	240	30	50
AS - 300 - 30 TF	S	300	30	50
AS - 350 - 30 TF	S	350	30	40
AS - 400 - 30 TF	S	400	30	40

SUCTION TUBES WITH UNIVERSAL CONNECTOR (T)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 T	S	150	25	100
AS - 180 - 25 T	S	180	25	100
AS - 200 - 25 T	S	200	25	100
AS - 210 - 25 T	S	210	25	100
AS - 240 - 25 T	S	240	25	75
AS - 300 - 25 T	S	300	25	60
AS - 350 - 25 T	S	350	25	50
AS - 400 - 25 T	S	400	25	50
AS - 150 - 30 T	S	150	30	50
AS - 180 - 30 T	S	180	30	50
AS - 200 - 30 T	S	200	30	50
AS - 210 - 30 T	S	210	30	50
AS - 240 - 30 T	S	240	30	50
AS - 300 - 30 T	S	300	30	50
AS - 350 - 30 T	S	350	30	40
AS - 400 - 30 T	S	400	30	40

SUCTION TUBES WITH FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 F	S	150	25	100
AS - 180 - 25 F	S	180	25	100
AS - 200 - 25 F	S	200	25	100
AS - 210 - 25 F	S	210	25	100
AS - 240 - 25 F	S	240	25	75
AS - 300 - 25 F	S	300	25	60
AS - 350 - 25 F	S	350	25	50
AS - 400 - 25 F	S	400	25	50
AS - 150 - 30 F	S	150	30	50
AS - 180 - 30 F	S	180	30	50
AS - 200 - 30 F	S	200	30	50
AS - 210 - 30 F	S	210	30	50
AS - 240 - 30 F	S	240	30	50
AS - 300 - 30 F	S	300	30	50
AS - 350 - 30 F	S	350	30	40
AS - 400 - 30 F	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITHOUT ATTACHMENTS




REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25	S	150	25	100
AS - 180 - 25	S	180	25	100
AS - 200 - 25	S	200	25	100
AS - 210 - 25	S	210	25	100
AS - 240 - 25	S	240	25	75
AS - 300 - 25	S	300	25	60
AS - 350 - 25	S	350	25	50
AS - 400 - 25	S	400	25	50
AS - 150 - 30	S	150	30	75
AS - 180 - 30	S	180	30	75
AS - 200 - 30	S	200	30	50
AS - 210 - 30	S	210	30	50
AS - 240 - 30	S	240	30	50
AS - 300 - 30	S	300	30	50
AS - 350 - 30	S	350	30	40
AS - 400 - 30	S	400	30	40

SUCTION TUBES WITH 2 UNIVERSAL CONNECTORS (TT)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TT	S	150	25	100
AS - 180 - 25 TT	S	180	25	100
AS - 200 - 25 TT	S	200	25	100
AS - 210 - 25 TT	S	210	25	100
AS - 240 - 25 TT	S	240	25	75
AS - 300 - 25 TT	S	300	25	60
AS - 350 - 25 TT	S	350	25	50
AS - 400 - 25 TT	S	400	25	50
AS - 150 - 30 TT	S	150	30	50
AS - 180 - 30 TT	S	180	30	50
AS - 200 - 30 TT	S	200	30	50
AS - 210 - 30 TT	S	210	30	50
AS - 240 - 30 TT	S	240	30	50
AS - 300 - 30 TT	S	300	30	50
AS - 350 - 30 TT	S	350	30	40
AS - 400 - 30 TT	S	400	30	40

SUCTION TUBES WITH UNIVERSAL CONNECTOR (T) AND LOOSE FINGERTIP (FL)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TFL	S	150	25	100
AS - 180 - 25 TFL	S	180	25	100
AS - 200 - 25 TFL	S	200	25	100
AS - 210 - 25 TFL	S	210	25	100
AS - 240 - 25 TFL	S	240	25	75
AS - 300 - 25 TFL	S	300	25	60
AS - 350 - 25 TFL	S	350	25	50
AS - 400 - 25 TFL	S	400	25	50
AS - 150 - 30 TFL	S	150	30	50
AS - 180 - 30 TFL	S	180	30	50
AS - 200 - 30 TFL	S	200	30	50
AS - 210 - 30 TFL	S	210	30	50
AS - 240 - 30 TFL	S	240	30	50
AS - 300 - 30 TFL	S	300	30	50
AS - 350 - 30 TFL	S	350	30	40
AS - 400 - 30 TFL	S	400	30	40

SUCTION TUBES

SUCTION TUBES WITH 2 STRAIGHT ATTACHMENTS (LL)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 LL	S	150	25	100
AS - 180 - 25 LL	S	180	25	100
AS - 200 - 25 LL	S	200	25	100
AS - 210 - 25 LL	S	210	25	100
AS - 240 - 25 LL	S	240	25	75
AS - 300 - 25 LL	S	300	25	60
AS - 350 - 25 LL	S	350	25	50
AS - 400 - 25 LL	S	400	25	50
AS - 150 - 30 LL	S	150	30	50
AS - 180 - 30 LL	S	180	30	50
AS - 200 - 30 LL	S	200	30	50
AS - 210 - 30 LL	S	210	30	50
AS - 240 - 30 LL	S	240	30	50
AS - 300 - 30 LL	S	300	30	50
AS - 350 - 30 LL	S	350	30	40
AS - 400 - 30 LL	S	400	30	40

SUCTION TUBES WITH ANTI-KINK CONICAL CONNECTOR (Tk) AND CATHETER CONNECTOR (Tb)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkTb	S	150	25	100
AS - 180 - 25 TkTb	S	180	25	100
AS - 200 - 25 TkTb	S	200	25	100
AS - 210 - 25 TkTb	S	210	25	100
AS - 240 - 25 TkTb	S	240	25	75
AS - 300 - 25 TkTb	S	300	25	60
AS - 350 - 25 TkTb	S	350	25	50
AS - 400 - 25 TkTb	S	400	25	50
AS - 150 - 30 TkTb	S	150	30	50
AS - 180 - 30 TkTb	S	180	30	50
AS - 200 - 30 TkTb	S	200	30	50
AS - 210 - 30 TkTb	S	210	30	50
AS - 240 - 30 TkTb	S	240	30	50
AS - 300 - 30 TkTb	S	300	30	50
AS - 350 - 30 TkTb	S	350	30	40
AS - 400 - 30 TkTb	S	400	30	40


SUCTION TUBES WITH 2 ANTI-KINK CONICAL CONNECTORS (TkTk)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkTk	S	150	25	100
AS - 180 - 25 TkTk	S	180	25	100
AS - 200 - 25 TkTk	S	200	25	100
AS - 210 - 25 TkTk	S	210	25	100
AS - 240 - 25 TkTk	S	240	25	75
AS - 300 - 25 TkTk	S	300	25	60
AS - 350 - 25 TkTk	S	350	25	50
AS - 400 - 25 TkTk	S	400	25	50
AS - 150 - 30 TkTk	S	150	30	50
AS - 180 - 30 TkTk	S	180	30	50
AS - 200 - 30 TkTk	S	200	30	50
AS - 210 - 30 TkTk	S	210	30	50
AS - 240 - 30 TkTk	S	240	30	50
AS - 300 - 30 TkTk	S	300	30	50
AS - 350 - 30 TkTk	S	350	30	40
AS - 400 - 30 TkTk	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH CATHETER CONNECTOR (Tb) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TbF	S	150	25	100
AS - 180 - 25 TbF	S	180	25	100
AS - 200 - 25 TbF	S	200	25	100
AS - 210 - 25 TbF	S	210	25	100
AS - 240 - 25 TbF	S	240	25	75
AS - 300 - 25 TbF	S	300	25	60
AS - 350 - 25 TbF	S	350	25	50
AS - 400 - 25 TbF	S	400	25	50
AS - 150 - 30 TbF	S	150	30	50
AS - 180 - 30 TbF	S	180	30	50
AS - 200 - 30 TbF	S	200	30	50
AS - 210 - 30 TbF	S	210	30	50
AS - 240 - 30 TbF	S	240	30	50
AS - 300 - 30 TbF	S	300	30	50
AS - 350 - 30 TbF	S	350	30	40
AS - 400 - 30 TbF	S	400	30	40

SUCTION TUBES WITH 2 CATHETER CONNECTORS (TbTb)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TbTb	S	150	25	100
AS - 180 - 25 TbTb	S	180	25	100
AS - 200 - 25 TbTb	S	200	25	100
AS - 210 - 25 TbTb	S	210	25	100
AS - 240 - 25 TbTb	S	240	25	75
AS - 300 - 25 TbTb	S	300	25	60
AS - 350 - 25 TbTb	S	350	25	50
AS - 400 - 25 TbTb	S	400	25	50
AS - 150 - 30 TbTb	S	150	30	50
AS - 180 - 30 TbTb	S	180	30	50
AS - 200 - 30 TbTb	S	200	30	50
AS - 210 - 30 TbTb	S	210	30	50
AS - 240 - 30 TbTb	S	240	30	50
AS - 300 - 30 TbTb	S	300	30	50
AS - 350 - 30 TbTb	S	350	30	40
AS - 400 - 30 TbTb	S	400	30	40

SUCTION TUBES WITH ANTI-KINK CONICAL CONNECTOR (Tk) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkF	S	150	25	100
AS - 180 - 25 TkF	S	180	25	100
AS - 200 - 25 TkF	S	200	25	100
AS - 210 - 25 TkF	S	210	25	100
AS - 240 - 25 TkF	S	240	25	75
AS - 300 - 25 TkF	S	300	25	60
AS - 350 - 25 TkF	S	350	25	50
AS - 400 - 25 TkF	S	400	25	50
AS - 150 - 30 TkF	S	150	30	50
AS - 180 - 30 TkF	S	180	30	50
AS - 200 - 30 TkF	S	200	30	50
AS - 210 - 30 TkF	S	210	30	50
AS - 240 - 30 TkF	S	240	30	50
AS - 300 - 30 TkF	S	300	30	50
AS - 350 - 30 TkF	S	350	30	40
AS - 400 - 30 TkF	S	400	30	40

SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgF	S	150	25	100
AS - 180 - 25 TgF	S	180	25	100
AS - 200 - 25 TgF	S	200	25	100
AS - 210 - 25 TgF	S	210	25	100
AS - 240 - 25 TgF	S	240	25	75
AS - 300 - 25 TgF	S	300	25	60
AS - 350 - 25 TgF	S	350	25	50
AS - 400 - 25 TgF	S	400	25	50
AS - 150 - 30 TgF	S	150	30	50
AS - 180 - 30 TgF	S	180	30	50
AS - 200 - 30 TgF	S	200	30	50
AS - 210 - 30 TgF	S	210	30	50
AS - 240 - 30 TgF	S	240	30	50
AS - 300 - 30 TgF	S	300	30	50
AS - 350 - 30 TgF	S	350	30	40
AS - 400 - 30 TgF	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND STRAIGHT ATTACHMENT (L)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgL	S	150	25	100
AS - 180 - 25 TgL	S	180	25	100
AS - 200 - 25 TgL	S	200	25	100
AS - 210 - 25 TgL	S	210	25	100
AS - 240 - 25 TgL	S	240	25	75
AS - 300 - 25 TgL	S	300	25	60
AS - 350 - 25 TgL	S	350	25	50
AS - 400 - 25 TgL	S	400	25	50
AS - 150 - 30 TgL	S	150	30	50
AS - 180 - 30 TgL	S	180	30	50
AS - 200 - 30 TgL	S	200	30	50
AS - 210 - 30 TgL	S	210	30	50
AS - 240 - 30 TgL	S	240	30	50
AS - 300 - 30 TgL	S	300	30	50
AS - 350 - 30 TgL	S	350	30	40
AS - 400 - 30 TgL	S	400	30	40


SUCTION TUBES WITH 2 STEPPED CONICAL CONNECTORS (TgTg)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgTg	S	150	25	100
AS - 180 - 25 TgTg	S	180	25	100
AS - 200 - 25 TgTg	S	200	25	100
AS - 210 - 25 TgTg	S	210	25	100
AS - 240 - 25 TgTg	S	240	25	75
AS - 300 - 25 TgTg	S	300	25	60
AS - 350 - 25 TgTg	S	350	25	50
AS - 400 - 25 TgTg	S	400	25	50
AS - 150 - 30 TgTg	S	150	30	50
AS - 180 - 30 TgTg	S	180	30	50
AS - 200 - 30 TgTg	S	200	30	50
AS - 210 - 30 TgTg	S	210	30	50
AS - 240 - 30 TgTg	S	240	30	50
AS - 300 - 30 TgTg	S	300	30	50
AS - 350 - 30 TgTg	S	350	30	40
AS - 400 - 30 TgTg	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND UNIVERSAL CONNECTOR (T)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgT	S	150	25	100
AS - 180 - 25 TgT	S	180	25	100
AS - 200 - 25 TgT	S	200	25	100
AS - 210 - 25 TgT	S	210	25	100
AS - 240 - 25 TgT	S	240	25	75
AS - 300 - 25 TgT	S	300	25	60
AS - 350 - 25 TgT	S	350	25	50
AS - 400 - 25 TgT	S	400	25	50
AS - 150 - 30 TgT	S	150	30	50
AS - 180 - 30 TgT	S	180	30	50
AS - 200 - 30 TgT	S	200	30	50
AS - 210 - 30 TgT	S	210	30	50
AS - 240 - 30 TgT	S	240	30	50
AS - 300 - 30 TgT	S	300	30	50
AS - 350 - 30 TgT	S	350	30	40
AS - 400 - 30 TgT	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 Tg	S	150	25	100
AS - 180 - 25 Tg	S	180	25	100
AS - 200 - 25 Tg	S	200	25	100
AS - 210 - 25 Tg	S	210	25	100
AS - 240 - 25 Tg	S	240	25	75
AS - 300 - 25 Tg	S	300	25	60
AS - 350 - 25 Tg	S	350	25	50
AS - 400 - 25 Tg	S	400	25	50
AS - 150 - 30 Tg	S	150	30	50
AS - 180 - 30 Tg	S	180	30	50
AS - 200 - 30 Tg	S	200	30	50
AS - 210 - 30 Tg	S	210	30	50
AS - 240 - 30 Tg	S	240	30	50
AS - 300 - 30 Tg	S	300	30	50
AS - 350 - 30 Tg	S	350	30	40
AS - 400 - 30 Tg	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND ANTI-KINK CONICAL CONNECTOR (Tk)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgTk	S	150	25	100
AS - 180 - 25 TgTk	S	180	25	100
AS - 200 - 25 TgTk	S	200	25	100
AS - 210 - 25 TgTk	S	210	25	100
AS - 240 - 25 TgTk	S	240	25	75
AS - 300 - 25 TgTk	S	300	25	60
AS - 350 - 25 TgTk	S	350	25	50
AS - 400 - 25 TgTk	S	400	25	50
AS - 150 - 30 TgTk	S	150	30	50
AS - 180 - 30 TgTk	S	180	30	50
AS - 200 - 30 TgTk	S	200	30	50
AS - 210 - 30 TgTk	S	210	30	50
AS - 240 - 30 TgTk	S	240	30	50
AS - 300 - 30 TgTk	S	300	30	50
AS - 350 - 30 TgTk	S	350	30	40
AS - 400 - 30 TgTk	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND CONUS (K)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgK	S	150	25	100
AS - 180 - 25 TgK	S	180	25	100
AS - 200 - 25 TgK	S	200	25	100
AS - 210 - 25 TgK	S	210	25	100
AS - 240 - 25 TgK	S	240	25	75
AS - 300 - 25 TgK	S	300	25	60
AS - 350 - 25 TgK	S	350	25	50
AS - 400 - 25 TgK	S	400	25	50
AS - 150 - 30 TgK	S	150	30	50
AS - 180 - 30 TgK	S	180	30	50
AS - 200 - 30 TgK	S	200	30	50
AS - 210 - 30 TgK	S	210	30	50
AS - 240 - 30 TgK	S	240	30	50
AS - 300 - 30 TgK	S	300	30	50
AS - 350 - 30 TgK	S	350	30	40
AS - 400 - 30 TgK	S	400	30	40

BUBBLE TUBES



REF/non-sterile	length (cm)	I.D.	PCS
BS 3 - 50000	50	3	1
BS 5 - 50000	50	5	1
BS 7 - 50000	50	7	1

FINGERTIP (F)



REF	sterile	PCS
700/00201	x	100
REF	non-sterile	PCS
700/00202	x	3.000

SCHLAUCHVERBINDER STERIL

Erhältlich in verschiedenen Durchmessern für einen Anschluss eines weiteren Absaug Schlauches. Die Größenangabe bezieht sich auf den Außendurchmesser an der kleinsten Stelle des Verbinders. Dabei können die beiden Seiten des Konnektors unterschiedliche Durchmesser aufweisen.



SCHLAUCHVERBINDER

REF 855/13107

Material: PA, weiß

Größe: 10/13

VPE: 250 Stk.



SCHLAUCHVERBINDER

REF 855/40470

Material: PA, weiß

Größe: 4/4

VPE: 250 Stk.



SCHLAUCHVERBINDER

REF 855/40770

Material: PA, weiß

Größe: 4/7

VPE: 250 Stk.



SCHLAUCHVERBINDER

REF 855/70770

Material: PA, weiß

Größe: 7/7

VPE: 250 Stk.



SCHLAUCHVERBINDER

REF 855/71070

Material: PA, weiß

Größe: 7/10

VPE: 250 Stk.



SCHLAUCHVERBINDER UNSTERIL



SCHLAUCHVERBINDER

REF 1010/13106

Material: PA, weiß

Größe: 10/13

VPE: 100 Stk.



SCHLAUCHVERBINDER

REF 1010/40460

Material: PA, weiß

Größe: 4/4

VPE: 100 Stk.



SCHLAUCHVERBINDER

REF 1010/40770

Material: PA, weiß

Größe: 4/7

VPE: 100 Stk.



SCHLAUCHVERBINDER

REF 1010/70760

Material: PA, weiß

Größe: 7/7

VPE: 100 Stk.



SCHLAUCHVERBINDER

REF 1010/71060

Material: PA, weiß

Größe: 7/10

VPE: 100 Stk.

INQUIRY ORDER

Company

Name

Address

Zip, City

Fon

Fax

E-Mail

REF	Product name	Number (cartons)

Date, signature



SEM  **Plastomed**



SEM-Plastomed GmbH
Medizin- und Kunststofftechnik
Schulstr. 6
D-57612 Obererbach
Fon: +49 (0) 2681 - 95 07 80
Fax: +49 (0) 2681 - 95 07 818
mail@plastomed.eu
www.plastomed.eu

