



Plastomed

Systems for your life...



SUCTION, BUBBLE & SILICON TUBES



Innovative plastic products and systems for medicine and technology

Product Description

As one of the leading European suppliers, we produce a wide-ranging & flexible assortment of medical tube systems. This extensive variety is made up through the application of distinctive attachment & adapter variations, tube types & lengths, tailor-made in accordance with your individual specifications.

All tube types are available individually wrapped in Peel Pack, either sterile or non-sterile.

Alternatively, we also pack your products singly into PE bags. They can then be integrated together and then package the products to form individual OP sets.

Suction hoses have to stay free from kinks in order to guarantee 100% perfect function. The right balance between flexibility & stability for our standard tubes is achieved through an optimised ratio of overall diameter to wall thickness. All tubes are also available as soft-tubes (soft material), just as in the case of standard tubes, also free from Latex and DEHP-free.

Our tube attachments / adapters guarantee safe connection / disconnection of tubes and are suitable for all suction-related units, bags and systems.

Notes when ordering

The following pages will give you an overview through our standard types. Furthermore, we manufacture & supply suction tubes in customised lengths, diameters and connectors- all in line with specialist client wishes.

Some examples of article composition

REF:	AS-180-25 TFS
AS:	Suction tubes
180:	Tube lengths in cm
25:	CH sizes for tubes
T:	Universal funnel
F:	Fingertip
S:	sterile

Application fields

ENT
Home-Care
Oral surgery
Orthopaedics
Surgery in general
Liposuction
Endoscopy
Gynaecology
Heart / Thorax surgery
Intensive care
Tracheostomy

Contents

Tube attachments/ Variations	page
Universal connectors	4
Stepped conical connectors	4
Anti-kink conical connectors	4
Catheter connectors	4
Straight attachments	5
Fingertip	5
Stepped connector / Conus	5
Suction tubes	
Suction tubes with universal connector and fingertip	6
Suction tubes with universal connector	6
Suction tubes with fingertip	6
Suction tubes without attachments	7
Suction tubes with 2 universal connectors	7
Suction tubes with universal connector and loose fingertip	7
Suction tubes with 2 straight attachments	8
Suction tubes with anti-kink conical connector and catheter connector	8
Suction tubes with 2 anti-kink conical connectors	8
Suction tubes with catheter connector and fingertip	9
Suction tubes with 2 catheter connectors	9
Suction tubes with anti-kink conical connector and fingertip	9
Suction tubes with stepped conical connector and fingertip	10
Suction tubes with stepped conical connector and straight attachment	10
Suction tubes with 2 stepped conical connectors	10
Suction tubes with stepped conical connector and universal connector	11
Suction tubes with stepped conical connector	11
Suction tubes with stepped conical connector and anti-kink conical connector	11
Suction tubes with stepped conical connector and conus	12
Bubble tubes	12
Fingertip	12
Silicon tubes and sets	
Silicon tubes	14
Ordering information	
Order formula	15

TUBE ATTACHMENTS / VARIATIONS



UNIVERSAL CONNECTORS (T)

Varied application onto suction pumps or bags with small connections. Tubes with universal funnel find use in fields such as in home-care, intensive care stations, OP theatres, tracheotomy or in ENT.



STEPPED CONICAL CONNECTORS (Tg)

Varied application onto suction pumps or bags with small and large connections. The anti-kink feature ensures a premature collapse of the tube is prevented under loading, because the adapter bends. Tubes with stepped conical connector find use in areas such as in OP theatres.



ANTI-KINK CONICAL CONNECTORS (Tk)


Varied application onto suction pumps or bags with small and large connections. The anti-kink feature ensures a premature collapse of the tube is prevented under loading, because the adapter bends. Tubes with anti-kink conical connector find use in areas such as in OP theatres.



CATHETER CONNECTORS (Tb)

Varied application onto suction pumps or bags with small and large connections. Tubes with catheter adapter find use in areas such as in in home-care, intensive care stations, OP theatres, tracheotomy or ENT.

TUBE ATTACHMENTS / VARIATIONS



STRAIGHT ATTACHMENTS (L)

Varied application onto suction pumps or bags with small and large connections. Tubes with straight attachments find use in areas such as in OP theatres.



FINGERTIP (F)

Interrupter attachment (fingertip) purpose: regulation of the suction flow of fluids and OP residues.



STEPPED CONNECTOR / CONUS (K)

The connector attaches the suction tube as extension to the suction unit with the thorax drainage. Tubes with stepped connectors find application in OP theatres.


SUCTION TUBES

SUCTION TUBES WITH UNIVERSAL CONNECTOR (T) AND FINGERTIP (F)




REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TF	S	150	25	100
AS - 180 - 25 TF	S	180	25	100
AS - 200 - 25 TF	S	200	25	100
AS - 210 - 25 TF	S	210	25	100
AS - 240 - 25 TF	S	240	25	75
AS - 300 - 25 TF	S	300	25	60
AS - 350 - 25 TF	S	350	25	50
AS - 400 - 25 TF	S	400	25	50
AS - 150 - 30 TF	S	150	30	50
AS - 180 - 30 TF	S	180	30	50
AS - 200 - 30 TF	S	200	30	50
AS - 210 - 30 TF	S	210	30	50
AS - 240 - 30 TF	S	240	30	50
AS - 300 - 30 TF	S	300	30	50
AS - 350 - 30 TF	S	350	30	40
AS - 400 - 30 TF	S	400	30	40

SUCTION TUBES WITH UNIVERSAL CONNECTOR (T)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 T	S	150	25	100
AS - 180 - 25 T	S	180	25	100
AS - 200 - 25 T	S	200	25	100
AS - 210 - 25 T	S	210	25	100
AS - 240 - 25 T	S	240	25	75
AS - 300 - 25 T	S	300	25	60
AS - 350 - 25 T	S	350	25	50
AS - 400 - 25 T	S	400	25	50
AS - 150 - 30 T	S	150	30	50
AS - 180 - 30 T	S	180	30	50
AS - 200 - 30 T	S	200	30	50
AS - 210 - 30 T	S	210	30	50
AS - 240 - 30 T	S	240	30	50
AS - 300 - 30 T	S	300	30	50
AS - 350 - 30 T	S	350	30	40
AS - 400 - 30 T	S	400	30	40


SUCTION TUBES WITH FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 F	S	150	25	100
AS - 180 - 25 F	S	180	25	100
AS - 200 - 25 F	S	200	25	100
AS - 210 - 25 F	S	210	25	100
AS - 240 - 25 F	S	240	25	75
AS - 300 - 25 F	S	300	25	60
AS - 350 - 25 F	S	350	25	50
AS - 400 - 25 F	S	400	25	50
AS - 150 - 30 F	S	150	30	50
AS - 180 - 30 F	S	180	30	50
AS - 200 - 30 F	S	200	30	50
AS - 210 - 30 F	S	210	30	50
AS - 240 - 30 F	S	240	30	50
AS - 300 - 30 F	S	300	30	50
AS - 350 - 30 F	S	350	30	40
AS - 400 - 30 F	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITHOUT ATTACHMENTS




REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25	S	150	25	100
AS - 180 - 25	S	180	25	100
AS - 200 - 25	S	200	25	100
AS - 210 - 25	S	210	25	100
AS - 240 - 25	S	240	25	75
AS - 300 - 25	S	300	25	60
AS - 350 - 25	S	350	25	50
AS - 400 - 25	S	400	25	50
AS - 150 - 30	S	150	30	75
AS - 180 - 30	S	180	30	75
AS - 200 - 30	S	200	30	50
AS - 210 - 30	S	210	30	50
AS - 240 - 30	S	240	30	50
AS - 300 - 30	S	300	30	50
AS - 350 - 30	S	350	30	40
AS - 400 - 30	S	400	30	40

SUCTION TUBES WITH 2 UNIVERSAL CONNECTORS (TT)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TT	S	150	25	100
AS - 180 - 25 TT	S	180	25	100
AS - 200 - 25 TT	S	200	25	100
AS - 210 - 25 TT	S	210	25	100
AS - 240 - 25 TT	S	240	25	75
AS - 300 - 25 TT	S	300	25	60
AS - 350 - 25 TT	S	350	25	50
AS - 400 - 25 TT	S	400	25	50
AS - 150 - 30 TT	S	150	30	50
AS - 180 - 30 TT	S	180	30	50
AS - 200 - 30 TT	S	200	30	50
AS - 210 - 30 TT	S	210	30	50
AS - 240 - 30 TT	S	240	30	50
AS - 300 - 30 TT	S	300	30	50
AS - 350 - 30 TT	S	350	30	40
AS - 400 - 30 TT	S	400	30	40


SUCTION TUBES WITH UNIVERSAL CONNECTOR (T) AND LOOSE FINGERTIP (FL)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TFL	S	150	25	100
AS - 180 - 25 TFL	S	180	25	100
AS - 200 - 25 TFL	S	200	25	100
AS - 210 - 25 TFL	S	210	25	100
AS - 240 - 25 TFL	S	240	25	75
AS - 300 - 25 TFL	S	300	25	60
AS - 350 - 25 TFL	S	350	25	50
AS - 400 - 25 TFL	S	400	25	50
AS - 150 - 30 TFL	S	150	30	50
AS - 180 - 30 TFL	S	180	30	50
AS - 200 - 30 TFL	S	200	30	50
AS - 210 - 30 TFL	S	210	30	50
AS - 240 - 30 TFL	S	240	30	50
AS - 300 - 30 TFL	S	300	30	50
AS - 350 - 30 TFL	S	350	30	40
AS - 400 - 30 TFL	S	400	30	40

SUCTION TUBES

SUCTION TUBES WITH 2 STRAIGHT ATTACHMENTS (LL)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 LL	S	150	25	100
AS - 180 - 25 LL	S	180	25	100
AS - 200 - 25 LL	S	200	25	100
AS - 210 - 25 LL	S	210	25	100
AS - 240 - 25 LL	S	240	25	75
AS - 300 - 25 LL	S	300	25	60
AS - 350 - 25 LL	S	350	25	50
AS - 400 - 25 LL	S	400	25	50
AS - 150 - 30 LL	S	150	30	50
AS - 180 - 30 LL	S	180	30	50
AS - 200 - 30 LL	S	200	30	50
AS - 210 - 30 LL	S	210	30	50
AS - 240 - 30 LL	S	240	30	50
AS - 300 - 30 LL	S	300	30	50
AS - 350 - 30 LL	S	350	30	40
AS - 400 - 30 LL	S	400	30	40

SUCTION TUBES WITH ANTI-KINK CONICAL CONNECTOR (Tk) AND CATHETER CONNECTOR (Tb)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkTb	S	150	25	100
AS - 180 - 25 TkTb	S	180	25	100
AS - 200 - 25 TkTb	S	200	25	100
AS - 210 - 25 TkTb	S	210	25	100
AS - 240 - 25 TkTb	S	240	25	75
AS - 300 - 25 TkTb	S	300	25	60
AS - 350 - 25 TkTb	S	350	25	50
AS - 400 - 25 TkTb	S	400	25	50
AS - 150 - 30 TkTb	S	150	30	50
AS - 180 - 30 TkTb	S	180	30	50
AS - 200 - 30 TkTb	S	200	30	50
AS - 210 - 30 TkTb	S	210	30	50
AS - 240 - 30 TkTb	S	240	30	50
AS - 300 - 30 TkTb	S	300	30	50
AS - 350 - 30 TkTb	S	350	30	40
AS - 400 - 30 TkTb	S	400	30	40


SUCTION TUBES WITH 2 ANTI-KINK CONICAL CONNECTORS (TkTk)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkTk	S	150	25	100
AS - 180 - 25 TkTk	S	180	25	100
AS - 200 - 25 TkTk	S	200	25	100
AS - 210 - 25 TkTk	S	210	25	100
AS - 240 - 25 TkTk	S	240	25	75
AS - 300 - 25 TkTk	S	300	25	60
AS - 350 - 25 TkTk	S	350	25	50
AS - 400 - 25 TkTk	S	400	25	50
AS - 150 - 30 TkTk	S	150	30	50
AS - 180 - 30 TkTk	S	180	30	50
AS - 200 - 30 TkTk	S	200	30	50
AS - 210 - 30 TkTk	S	210	30	50
AS - 240 - 30 TkTk	S	240	30	50
AS - 300 - 30 TkTk	S	300	30	50
AS - 350 - 30 TkTk	S	350	30	40
AS - 400 - 30 TkTk	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH CATHETER CONNECTOR (Tb) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TbF	S	150	25	100
AS - 180 - 25 TbF	S	180	25	100
AS - 200 - 25 TbF	S	200	25	100
AS - 210 - 25 TbF	S	210	25	100
AS - 240 - 25 TbF	S	240	25	75
AS - 300 - 25 TbF	S	300	25	60
AS - 350 - 25 TbF	S	350	25	50
AS - 400 - 25 TbF	S	400	25	50
AS - 150 - 30 TbF	S	150	30	50
AS - 180 - 30 TbF	S	180	30	50
AS - 200 - 30 TbF	S	200	30	50
AS - 210 - 30 TbF	S	210	30	50
AS - 240 - 30 TbF	S	240	30	50
AS - 300 - 30 TbF	S	300	30	50
AS - 350 - 30 TbF	S	350	30	40
AS - 400 - 30 TbF	S	400	30	40

SUCTION TUBES WITH 2 CATHETER CONNECTORS (TbTb)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TbTb	S	150	25	100
AS - 180 - 25 TbTb	S	180	25	100
AS - 200 - 25 TbTb	S	200	25	100
AS - 210 - 25 TbTb	S	210	25	100
AS - 240 - 25 TbTb	S	240	25	75
AS - 300 - 25 TbTb	S	300	25	60
AS - 350 - 25 TbTb	S	350	25	50
AS - 400 - 25 TbTb	S	400	25	50
AS - 150 - 30 TbTb	S	150	30	50
AS - 180 - 30 TbTb	S	180	30	50
AS - 200 - 30 TbTb	S	200	30	50
AS - 210 - 30 TbTb	S	210	30	50
AS - 240 - 30 TbTb	S	240	30	50
AS - 300 - 30 TbTb	S	300	30	50
AS - 350 - 30 TbTb	S	350	30	40
AS - 400 - 30 TbTb	S	400	30	40


SUCTION TUBES WITH ANTI-KINK CONICAL CONNECTOR (Tk) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TkF	S	150	25	100
AS - 180 - 25 TkF	S	180	25	100
AS - 200 - 25 TkF	S	200	25	100
AS - 210 - 25 TkF	S	210	25	100
AS - 240 - 25 TkF	S	240	25	75
AS - 300 - 25 TkF	S	300	25	60
AS - 350 - 25 TkF	S	350	25	50
AS - 400 - 25 TkF	S	400	25	50
AS - 150 - 30 TkF	S	150	30	50
AS - 180 - 30 TkF	S	180	30	50
AS - 200 - 30 TkF	S	200	30	50
AS - 210 - 30 TkF	S	210	30	50
AS - 240 - 30 TkF	S	240	30	50
AS - 300 - 30 TkF	S	300	30	50
AS - 350 - 30 TkF	S	350	30	40
AS - 400 - 30 TkF	S	400	30	40

SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND FINGERTIP (F)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgF	S	150	25	100
AS - 180 - 25 TgF	S	180	25	100
AS - 200 - 25 TgF	S	200	25	100
AS - 210 - 25 TgF	S	210	25	100
AS - 240 - 25 TgF	S	240	25	75
AS - 300 - 25 TgF	S	300	25	60
AS - 350 - 25 TgF	S	350	25	50
AS - 400 - 25 TgF	S	400	25	50
AS - 150 - 30 TgF	S	150	30	50
AS - 180 - 30 TgF	S	180	30	50
AS - 200 - 30 TgF	S	200	30	50
AS - 210 - 30 TgF	S	210	30	50
AS - 240 - 30 TgF	S	240	30	50
AS - 300 - 30 TgF	S	300	30	50
AS - 350 - 30 TgF	S	350	30	40
AS - 400 - 30 TgF	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND STRAIGHT ATTACHMENT (L)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgL	S	150	25	100
AS - 180 - 25 TgL	S	180	25	100
AS - 200 - 25 TgL	S	200	25	100
AS - 210 - 25 TgL	S	210	25	100
AS - 240 - 25 TgL	S	240	25	75
AS - 300 - 25 TgL	S	300	25	60
AS - 350 - 25 TgL	S	350	25	50
AS - 400 - 25 TgL	S	400	25	50
AS - 150 - 30 TgL	S	150	30	50
AS - 180 - 30 TgL	S	180	30	50
AS - 200 - 30 TgL	S	200	30	50
AS - 210 - 30 TgL	S	210	30	50
AS - 240 - 30 TgL	S	240	30	50
AS - 300 - 30 TgL	S	300	30	50
AS - 350 - 30 TgL	S	350	30	40
AS - 400 - 30 TgL	S	400	30	40


SUCTION TUBES WITH 2 STEPPED CONICAL CONNECTORS (TgTg)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgTg	S	150	25	100
AS - 180 - 25 TgTg	S	180	25	100
AS - 200 - 25 TgTg	S	200	25	100
AS - 210 - 25 TgTg	S	210	25	100
AS - 240 - 25 TgTg	S	240	25	75
AS - 300 - 25 TgTg	S	300	25	60
AS - 350 - 25 TgTg	S	350	25	50
AS - 400 - 25 TgTg	S	400	25	50
AS - 150 - 30 TgTg	S	150	30	50
AS - 180 - 30 TgTg	S	180	30	50
AS - 200 - 30 TgTg	S	200	30	50
AS - 210 - 30 TgTg	S	210	30	50
AS - 240 - 30 TgTg	S	240	30	50
AS - 300 - 30 TgTg	S	300	30	50
AS - 350 - 30 TgTg	S	350	30	40
AS - 400 - 30 TgTg	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND UNIVERSAL CONNECTOR (T)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgT	S	150	25	100
AS - 180 - 25 TgT	S	180	25	100
AS - 200 - 25 TgT	S	200	25	100
AS - 210 - 25 TgT	S	210	25	100
AS - 240 - 25 TgT	S	240	25	75
AS - 300 - 25 TgT	S	300	25	60
AS - 350 - 25 TgT	S	350	25	50
AS - 400 - 25 TgT	S	400	25	50
AS - 150 - 30 TgT	S	150	30	50
AS - 180 - 30 TgT	S	180	30	50
AS - 200 - 30 TgT	S	200	30	50
AS - 210 - 30 TgT	S	210	30	50
AS - 240 - 30 TgT	S	240	30	50
AS - 300 - 30 TgT	S	300	30	50
AS - 350 - 30 TgT	S	350	30	40
AS - 400 - 30 TgT	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 Tg	S	150	25	100
AS - 180 - 25 Tg	S	180	25	100
AS - 200 - 25 Tg	S	200	25	100
AS - 210 - 25 Tg	S	210	25	100
AS - 240 - 25 Tg	S	240	25	75
AS - 300 - 25 Tg	S	300	25	60
AS - 350 - 25 Tg	S	350	25	50
AS - 400 - 25 Tg	S	400	25	50
AS - 150 - 30 Tg	S	150	30	50
AS - 180 - 30 Tg	S	180	30	50
AS - 200 - 30 Tg	S	200	30	50
AS - 210 - 30 Tg	S	210	30	50
AS - 240 - 30 Tg	S	240	30	50
AS - 300 - 30 Tg	S	300	30	50
AS - 350 - 30 Tg	S	350	30	40
AS - 400 - 30 Tg	S	400	30	40

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND ANTI-KINK CONICAL CONNECTOR (Tk)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgTk	S	150	25	100
AS - 180 - 25 TgTk	S	180	25	100
AS - 200 - 25 TgTk	S	200	25	100
AS - 210 - 25 TgTk	S	210	25	100
AS - 240 - 25 TgTk	S	240	25	75
AS - 300 - 25 TgTk	S	300	25	60
AS - 350 - 25 TgTk	S	350	25	50
AS - 400 - 25 TgTk	S	400	25	50
AS - 150 - 30 TgTk	S	150	30	50
AS - 180 - 30 TgTk	S	180	30	50
AS - 200 - 30 TgTk	S	200	30	50
AS - 210 - 30 TgTk	S	210	30	50
AS - 240 - 30 TgTk	S	240	30	50
AS - 300 - 30 TgTk	S	300	30	50
AS - 350 - 30 TgTk	S	350	30	40
AS - 400 - 30 TgTk	S	400	30	40


SUCTION TUBES

SUCTION TUBES WITH STEPPED CONICAL CONNECTOR (Tg) AND CONUS (K)



REF/non-sterile	sterile	length (cm)	CH-size	PCS
AS - 150 - 25 TgK	S	150	25	100
AS - 180 - 25 TgK	S	180	25	100
AS - 200 - 25 TgK	S	200	25	100
AS - 210 - 25 TgK	S	210	25	100
AS - 240 - 25 TgK	S	240	25	75
AS - 300 - 25 TgK	S	300	25	60
AS - 350 - 25 TgK	S	350	25	50
AS - 400 - 25 TgK	S	400	25	50
AS - 150 - 30 TgK	S	150	30	50
AS - 180 - 30 TgK	S	180	30	50
AS - 200 - 30 TgK	S	200	30	50
AS - 210 - 30 TgK	S	210	30	50
AS - 240 - 30 TgK	S	240	30	50
AS - 300 - 30 TgK	S	300	30	50
AS - 350 - 30 TgK	S	350	30	40
AS - 400 - 30 TgK	S	400	30	40

BUBBLE TUBES



REF/non-sterile	length (cm)	I.D.	PCS
BS 3 - 50000	50	3	1
BS 5 - 50000	50	5	1
BS 7 - 50000	50	7	1

FINGERTIP (F)



REF	sterile	PCS
700/00201	x	100

SILICON TUBES AND SETS



SILICON TUBES AND SETS

SILICON TUBES

Silicon tubes, non-sterile

I.D. (mm)	O.D. (mm)	I.D. (inch)	O.D. (inch)	REF	length (cm)
1.0	2.0	.039"	.079"	SI 105	50
1.0	3.0	.039"	.118"	SI 110	50
2.0	4.0	.079"	.157"	SI 210	50
3.0	5.0	.118"	.197"	SI 310	50
3.0	6.0	.118"	.236"	SI 315	50
4.0	6.0	.157"	.236"	SI 410	50
4.0	7.0	.157"	.276"	SI 415	50
4.0	8.0	.157"	.315"	SI 420	50
5.0	7.0	.197"	.276"	SI 510	50
5.0	8.0	.197"	.315"	SI 515	50
5.0	9.0	.197"	.354"	SI 520	50
6.0	8.0	.236"	.315"	SI 610	50
6.0	9.0	.236"	.354"	SI 615	50
6.0	10.0	.236"	.394"	SI 620	50
6.0	12.0	.236"	.472"	SI 630	50
7.0	10.0	.276"	.394"	SI 715	50
7.0	11.0	.276"	.433"	SI 720	50
7.0	13.0	.276"	.512"	SI 730	50
8.0	12.0	.315"	.472"	SI 820	50
8.0	14.0	.315"	.551"	SI 830	50
10.0	13.0	.315"	.472"	SI 820	50

Optional extras include silicon tubes, cut into fixed lengths, confectioned with fingertip or filter, packaged, labelled and sterilised with ETO gas if so desired.



INQUIRY ☐ ORDER

Company

Name

Address

Zip, City


Fon

Fax

E-Mail

REF	Product name	Number (cartons)

Date, signature



SEM - Plastomed GmbH
Schulstr. 6
D-57612 Obererbach
Fon: +49 (0) 2681 - 95 07 80
Fax: +49 (0) 2681 - 95 07 818
mail@plastomed.eu
www.plastomed.eu

