



Suction, bladder and silicone tubings

Intended purpose

Suction, bladder and silicone tubings are used as connecting tubes between suction device, suction container and suction tips or suction catheter. Thus the tubes are transporting (sucking) blood, secrets, fats, tumescence solutions, etc. from the patient to the device.

This usage can be sustained only, when the tubes are installed without kinks.

These products are **non-invasive single use medical devices class I** (sterile, non-sterile) for transient usage (Normally intended for continuous use for less than 60 minutes) in connection with an **active** medical device. For hygienic reasons the device is to be exchanged after usage.

Classification

According to rule 1:

All non-invasive devices are in Class I, unless one of the rules set out hereinafter applies. **(93/42 EWG, Appendix IX, III Classification)**

Reason for class 1:

The tubes are used to transport (suck) blood, secrets, fats, tumescence solutions, etc. from the patient to the device.