



Plastomed

Systeme für Ihr Leben...

Surgical suction sets / Handpieces

Intended purpose

Surgical suction sets

The suction cannula / handpieces (Yankauer) are used as a suction tip. The surgical suction sets, consisting of handpiece and suction tube, are used as a connecting tube between the suction device, the cup and the patient. The surgical suction sets have the function of operating residues, blood, secretions, etc. from the patient to the device to pass through (suck off). The operation can only be maintained if the hoses are installed kink-free, to ensure the suction.

Yankauer / handpieces

The Yankauer / handpieces are for suction of blood, surgery residues and secretions.

Classification

According to rule 6:

All Temporary Surgery invasive Products are in Class IIa. **(93/42 EWG, Appendix IX, III Classification)**

Reason for class IIa:

The product is used to connect to an **active class IIa or higher medical device** in surgically invasive procedures and is used to **lead through blood or other body fluids.**